

Library and templates

MolDraw Docs · LabCanvas · 1 pages · moldraw.com/docs/

Library and templates

Short answer: Open Library in the ring bar to place ready-made structures. It has functional groups, ligands, amino acids, polymers, COFs, MOFs and named reactions.

- 1 Click **Library** at the end of the ring bar.
- 2 Search or pick a category.
- 3 Click an item, then click the canvas to place it.

Library is not the same as **My designs**. My designs holds your own saved drawings.

Questions and answers

Where are amino acid templates?

Source: <https://www.moldraw.com/docs/library/how-to-use-the-structure-library-and-templates/>